Mode d'emploi
Instruction for use
Verarbeitungsanleitung
Istruzioni d'uso
Modo de empleo

PX EXTRAVEST

PX EXTRAVEST: Mode d'emploi

Description

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PX EXTRAVEST: Instruction for use

Description

PX Extravest is a phosphate-bonded, graphite-free speed cast investment with exceptional properties. It is suitable for casting all types of dental alloy and can be preheated up to 1100°C. It can be used for crown and bridge work using precious and non-precious alloys, metal denture bases, fixed/removable restorations and injection-moulded porcelains (e.g. Empress®, Ivoclar).

Technical data

-Mixing liquid:	Expansion liquid (1200 ml)
-Mixing ratio powder/liquid :	100 g/25 ml (Crown and bridge), 160 g/40 ml. 100 g/23 ml (CrCo), 160 g/37 ml.
-Mixing under vacuum:	60 sec.
-Working time:	5–7 min.
-Working temperature:	17°C - 22°C
-Initial setting time complies with DIN 13919-B:	9 – 11 min.
-Setting expansion complies with DIN 13919-B:	appr. 2.00% (100% undiluted liquid)
-Recommended storage:	8°C – 12°C in a dry place

Crown and bridge work

Mix PX Extravest according to the instructions and pour it up to the patterns without vibration. Then vibrate gently into the patterns using a small instrument, if necessary. Fill the ring completely without any further vibration. Adhere strictly to the time given in the instructions for the speed technique.

Pouring models (CrCo)

Smear Vaseline on the inside of the overlapping section of the duplicating adhesive tape. Spray the duplicating mould with a silicone, then pour PX Extravest investment into the mould using light vibration. The model can be removed from the mould after approx. 40 minutes. Before waxing up, apply PX Extravest wax adhesive to areas of the model where preformed wax patterns are to be placed.

Speed casting

Preheat the furnace to 850°C and, exactly 15 minutes after mixing the investment, place the mould in the furnace with the sprue hole facing downwards. Adhere strictly to the prescribed time, which is measured from the first contact of the powder and liquid. Moulds should not be placed directly onto the floor of the furnace. We recommend using either a grooved tray for collecting wax or special stilts from a pottery craftshop. Caution! During the first 15 minutes the furnace fan and extractor should remain switched off and the furnace door should not be opened because of the risk of instant combustion. If a preheating temperature other than 850°C is required, the temperature can be adjusted 15 minutes after placing the mould in the furnace. Moulds are ready to cast after heat soaking for 45 – 60 min. depending on the size of the ring.

Conventional preheating

Allow the mould to set for at least 60 minutes, then place it in a cold furnace and heat to the required temperature at a heat rate of 5°C - 7°C per minute. Holding stages are not required when using PX Extravest.

Injection-moulded porcelains (e.g. Empress®)

Follow normal procedure and select a conventional linear preheating curve without holding stages. In many cases it is also possible to use the speed technique. The following recommendations for liquid concentrations for injection-moulded porcelains have been tried and tested in practice: 50% for inlays with 1 or 2 surfaces, 60% for inlays with 3 surfaces, 70% for the all-porcelain build-up technique and 80% for the all-porcelain staining technique. Carry out some test injections based on the above concentrations to determine the correct ratios for the injection-moulded porcelain used in your laboratory.

Please note

Technical recommendations are based on tests and findings from work in our development laboratory and can only be regarded as guidelines. PX DENTAL's products are subjected to strict quality controls. We reserve the right to make technical changes. Further information can be obtained from PX DENTAL technical services in La Chaux-de-Fonds, Switzerland (Phone: +41 32 924 21 20) or from our sales representatives.

Caution

Inhaling silica dust is a health risk. A suitable face mask should be worn.

IMPORTANT DATA

Mixing ratio Powder / Liquid

Crown and bridge work: 100g / 25 ml			Pouring models (CrCo): 100g / 23 ml	
ring size	powder	liquid	powder	liquid
1X	60 g	15 ml	60 g	14 ml
	100 g	25 ml	100 g	23 ml
3X	160 g	40 ml	160 g	37 ml
6X	320 g	80 ml	320 g	74 ml
9X	450 g	112 ml	450 g	104 ml

Concentration of the expansion liquid

(Max. expansion with undiluted liquid)

	Mixing ratio Powder / Liquid					
	Concentra- tion of liquid	1X = 60 g	100 g	3X = 160 g	6X = 320 g	9X = 450 g
High gold ceramic alloys	70-80%	11 ml / 4 ml	19 ml / 6 ml	30 ml / 10 ml	60 ml / 20 ml	84 ml / 28 ml
Gold reduced ceramic alloys	75-85%	12 ml / 3 ml	20 ml / 5 ml	32 ml / 8 ml	64 ml / 16 ml	90 ml / 22 ml
High gold casting alloys	60-65%	9 ml / 6 ml	16 ml / 9 ml	25 ml / 15 ml	50 ml / 30 ml	281 ml / 169 ml
Gold reduced casting alloys	65-70%	10 ml / 5 ml	17 ml / 8 ml	27 ml / 13 ml	54 ml / 26 ml	76 ml / 36 ml
CrCo/Fixed/removable appliances	75-85%	11 ml / 3 ml	18 ml / 5 ml	30 ml / 7 ml	60 ml / 14 ml	83 ml / 21 ml
Non-precious alloys	>90%	>13 ml / 2 ml	>22 ml / 3 ml	>36 ml / 4 ml	>72 ml / 8 ml	>100 ml / 12 ml
Injection-moulded porcelains	approx. 80%	12 ml / 3 ml	20 ml / 5 ml	32 ml / 8 ml	64 ml / 16 ml	90 ml / 22 ml

Steps to follow for investing

	Waiting time after mixing	Preheating of the furnace	Heating rate	Holding time at temperature	Remarks
Speed casting (max. 1100°C)	15 min! (important)	850°C	---	45-60 min.	During the first 15 minutes the furnace fan and extractor should remain switched off and the furnace door should not be opened because of the risk of instant combustion. If a preheating temperature other than 850°C is required, the temperature can be adjusted 15 minutes after placing the mould in the furnace.

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