PX N61**PXDENTAL SA**

CH-2074 Marin

Blanc Alliage dentaire céramo-metallique à base de Ni-Cr,
 Weiss Dental-Aufbrennlegierung auf Ni-Cr-Basis
 White Dental bonding alloy base Ni-Cr
 Blanco Aleación ceramo-metalica a base de Ni-Cr

% = Ni 61.0; Cr 26.0 ; Mo 11, Si 1.5 : <1% Fe, Ce, Al, Co

8.2 g/cm³

1280 – 1350°C

160 – 180 HV1

14.0 / 14.2 µm/m°C

ISO 22674 & 9693

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ISO 22674 & 9693

Instructions for Use

Alloy extra hard, for Crown & bridge technique and for cast partial dental frames. Free of beryllium and cadmium. Should not be use in patients suffering from nickel allergy.

Waxing up, spruing, investing

Follow the general prescriptions for bonding alloys.
The ultimate thickness of the cast coping should not be less than 0.3 mm.
Use only a high-heat phosphate-bonded investment material.

Preheating

Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
Keep the cylinder at **950° C for a soaking period 30 – 60 min** depending on ring size.

Crucibles and casting buttons

Use a ceramic crucible
Casting buttons through the melting of the alloy, part of the elements are lost. **Only a new material should be used.**

Melting

Induction casting machine: Allow for 1–2 seconds after the oxide layer has cracked open, then release the centrifuge.

Melting with an open flame (torch): Only torches with „multi orifice heads“ must be used, meaning that the heat must cover the largest possible area of the alloy being melted. Preheat the crucible in the furnace without alloy. When melting with a torch an oxide layer builds up over the coalesced ingots. Continue heating, until the alloy can be visibly moved by the flame under the oxide layer. **Release centrifuge before oxide layer breaks open.**
Bench cool the casting ring.

Treatment of the casting

Divesting: by blasting (grain size of 200-500 µm).
Finishing: use abrasive grinding points of your choice. Grind with ceramic bonded abrasive stones, tungsten-carbide and/or sintered diamond burs. Clean the frameworks in an ultrasonic unit or with a steam cleaner.

Oxidation




With **PX N61** an oxidation process is not required. If the manufacturer of the ceramic specifies, the oxide should be removed. Blasting (grain size of 125-250 µm) and the surface should be cleaned again.

Porcelain application

Porcelain applications are subject to the manufacturers instructions.
PX N61 is compatible with most known porcelain brands (high fusing).

Laser welding:

PX N61 is a carbon free alloy, very suitable for laser welding, due to their very low carbon contents

FR	DE	IT	EN	ES
 Consulter le mode d'emploi www.px dental.com	 Gebrauchsanweisung beachten www.px dental.com	 Consultare le istruzioni per l'uso www.px dental.com	 Consult instructions for use www.px dental.com	 Observe el modo de empleo www.px dental.com
REF	Katalognummer	Numero di catalogo	Catalogue number	Número de catálogo
LOT	Chargencode	Numero di lotto	Batch code	Número de lote
QTY	Quantität	Quantità	Quantity	Cantidad
				
				

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