Nacera[®] Zirconia

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INSTRUCTIONS FOR USE

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EN | INSTRUCTIONS FOR USE

EXPLANATION OF THE SYMBOLS USED

Symbol	Explanation
	Manufacturer
	Item number
	Batch
	Expiry date
	Observe the instructions for use
	Temperature limit
	Store in a dry place
	Prescription only
	Shrinkage factor
	Medical Device

In these instructions for use, the following products are referred to collectively as **Nacera® Blanks**:

- **Nacera® Pearl 1**
- **Nacera® Shell**
- **Nacera® Shell Technical**
- **Nacera® Pearl Shaded 16+2**
- **Nacera® Pearl Multi-Shade**
- **Nacera® Pearl Q³ Multi-Shade**
- **Nacera® Pearl Natural**

The physical properties of the individual products can be found in the technical data in the download area on our website:

www.nacera-medical.com/en

GENERAL HANDLING NOTES

Please read the instructions for use carefully before removing the **Nacera® Blank** from the packaging. They contain important information which is essential to ensure both the correct processing and the safety of patients and users.

Nacera® Blanks are manufactured and tested according to highest quality standards. In order to guaranty this level of quality during subsequent processing, the procedures described in the following must absolutely be adhered to.

Nacera® Blanks are supplied in a pre-sintered condition. They have limited stability or strength

and a residual porosity. Therefore, these blanks must be handled with care.

The blanks must be stored in their original packaging, at temperatures between 10 °C and 50 °C. They must not be subjected to impact or vibration. Contamination must absolutely be avoided. Make sure that the mounts are only handled with dry and clean hands or gloves and that they are under no circumstances contaminated with liquids (such as adhesives or marker pens).

If tooth stumps and metal mounts are heavily discoloured, we recommend our opaque material (**Nacera® Shell**).

WARNING

When using this product, the dentist must consider possible interaction between this medical product and other products already in place in the patient's mouth.

DANGER

The processing of blanks and finally-sintered mounts creates dust which might damage the lungs as well as irritating the eyes and skin. Therefore, processing may only be done

provided that the extractor system works correctly as well as wearing safety goggles and a filtering face piece. Please forward all of the information above to the dentist if you process this product for use in custom-made products within the framework of Directive 93/42 EC or Regulation 2017/245 (MDR).

Please also note the safety data sheets.

PURPOSE

Nacera® Blanks are zirconium oxide blanks for the production of fixed dentures, such as single crowns and bridges with up to 16* units (except **Nacera® Pearl Q³ Multi-Shade**). In the posterior tooth region, the span between the pillars must not exceed two units. A cantilever bridge of premolar size is permissible. The application is carried out by dental specialists.

*except Canada (up to 6 units)

Nacera® Pearl Q³ Multi-Shade blanks are Zirconium oxide blanks for the production of fixed dentures, such as single crowns and bridges with up to 3 units. In the posterior tooth region, the span between the pillars must not exceed a premolar width. The application is carried out by dental specialists.

INDICATION

Zirconium dioxide for the replacement of lost tooth substance (see indication table on page 58).

CONTRAINDICATIONS

For bruxism, the vertical dimension must be observed. In these cases, please consult with the dentist. In the event that there is an inadequate vertical dimension and that the preparation is unsuitable for an all-ceramic restoration, an alternative material must be chosen. Inlay bridges, enossal implants and root posts are other contraindications.

PREPARATION

Recommendations for the preparation include a marked chamfer or a rounded step.

- Minimum cutting depth at the preparation border 1 mm
- 1.5 to 2 mm of substance removed occlusal/incisal
- Edge radius 0.7 mm
- Preparation angle 6° to 8°

For bridge constructions, observe the parallelism and avoid negative steps. In general, please

observe the notes in the specialist literature.

NOTES ON CONSTRUCTION

Please place the construction in the layers of the **Nacera® Pearl Multi-Shade / Pearl Natural / Pearl Q³ Multi-Shade** blanks in so that the required colour / translucency is achieved.

Wall thickness **Nacera® Blanks**: For all-ceramic dentures, the wall thickness of the crowns in the sintered state must not be less than the corresponding values of table*. The minimum wall thickness must always be observed even after any adjustment of the mount by the dentist. The geometry of the connectors is of the greatest importance for the crack resistance of mounts made of zirconium oxide. Therefore, the cross-section of the connectors must be as large as possible.

For static reasons, the height of the connector is especially important. In order to avoid so-called „chipping“, the mounts should correspond to a reduced, fully anatomic shape of the crowns and bridges to ensure that the applied ceramic has the best possible support.

Please also refer to the processing instructions of the respective layering ceramic manufacturer.

The frameworks must be designed according to the general guidelines of digital dental technology.

*Wall thickness	
individual copings	bridges
0,4 mm	0,5 mm

Edge thickness	
individual copings	bridges
0,2 mm	0,2 mm

Mount dimensions		
	Nacera® blanks*	Nacera® Pearl Q ³ Multi-Shade
Intermediate segments	2	1 (premolar width)
Connector cross section	9 mm ² (posterior tooth region) 6 mm ² (anterior tooth region)	min. 9 mm ²
Cantilever bridge	1 (premolar width)	-
Connector cross section of cantilever bridge	12 mm ²	-

*Except **Nacera® Pearl Q³ Multi-Shade**

COLOURING THE BLANKS

Staining of the restorations (except **Nacera® Pearl Q³ Multi-Shade / Shell**) can be done with suitable staining liquids. Please observe the manufacturer's instructions! Drying of the restorations treated with Liquids, under red light or in the furnace according to the manufacturer's instructions, is absolutely necessary before final sintering.

MILLING STRATEGY

In your CAM software, please select „Common strategy“ for **Nacera® Shell** and for other **Nacera® Blanks**, please select the “highly translucent strategy“.

PROCESSING

Only those machines and tools which are approved for processing pre-sintered blanks made of zirconium oxide may be used for manufacturing dentures made of **Nacera® Blanks**.

Caution! For **Nacera® Pearl Multi-Shade / Pearl Q³ Multi-Shade / Pearl Natural** blanks, the printed side is the occlusal side.

If possible, neither coolant nor compressed

air should be used during processing. Please observe the instructions for use of your milling machine as well as the parameters of the CAD/CAM software.

Please separate the finished milled parts from the blanks carefully using suitable tools.


VISUAL INSPECTION

Before the milled mounts are processed further, they must be inspected for the following faults:


- Shiny areas on the surface (indicating a worn milling cutter)
- Material spalling (due to the milling strategy and milling cutter)
- Faulty mounts must not be processed further.

FINAL SINTERING

Sinter the **Nacera® Blanks** (except **Nacera® Pearl Q³ Multi-Shade**) in sintering program 1 (1500° C, see graphic sintering program 1). Sinter the **Nacera® Pearl Q³ Multi-Shade** blanks in sintering program 2 (1450° C, see graphic sintering program 2).

Sintering program 1

Sintering program 2

FURTHER PROCESSING OF THE FINALLY-SINTERED RESTORATIONS

If possible, the finally-sintered units should not be processed mechanically. However, if this is unavoidable, only water-cooled diamond-fitted tools in correct working order may be used. The interdental connection points should not be polished. In principle, the basal polishing of these connection points must always be avoided for stability reasons. Also in the case of abutments, sharp edges should be avoided and instead rounding is to be aimed for.

VENEERING

Veneering is done with commonly available veneer ceramics that are approved for zirconium oxide according to manufacturer specifications.

INDIVIDUALISATION

For individualising the restoration, painting techniques as well as cut-back and layering techniques, or a combination of both, are suitable.

GRINDING

To protect the antagonists and for material engineering reasons, ground contact points and surfaces must either be high-gloss polished after the try-in and/or glaze-fired.

INCORPORATING

State-of-the-art dental technology prefers self-adhesive and adhesive fastening over conventional fastening.

The fastening logs of the cement manufacturers must be observed.

DISPOSAL CONSIDERATIONS

Disposal of the product must be in accordance

with the regulatory requirements. Do not allow it to be released into the sewer system or the environment. Do not dispose of it together with the household waste.

REPORTING OF INCIDENTS

If serious incidents occur related to the product, report them to the manufacturer and your relevant regional authority.

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